



Leihu Medical Product Catalog

Product		Medical external fixation device				
Company		Shenzhen Leihu Medical Equipment Manufacturing Co., LTD.				
Position	Applications	Model	Before shaping	After shaping	Size selection reference	Size reference diagram
Hand	Wrist joint Mid-to-distal ulna/radius	LH-A101-S (Small)			A=22CM B=20CM C=16CM	
		LH-A101-M (Medium)			A=27CM B=25CM C=21CM	
		LH-A101-L (Large)			A=32CM B=30CM C=25CM	
		LH-A101-XL (Extra Large)			A=35CM B=32CM C=27CM	
Elbow	Mid-to-proximal ulna/radius	LH-A102-S (Small)			A=27CM B=15.5CM C=12CM	
		LH-A102-M (Medium)			A=31CM B=16.9CM C=15CM	
	Elbow joint	LH-A102-L (Large)			A=36CM B=19.2CM C=17CM	
	Distal humerus	LH-A102-XL (Extra Large)			A=40CM B=20.5CM C=19CM	
Foot108°	Ankle joint (108°)	LH-A201-S (Small)			A=14CM B=19CM C=20CM	
		LH-A201-M (Medium)			A=17CM B=24CM C=25CM	
	Distal tibia/fibula	LH-A201-L (Large)			A=19CM B=32CM C=28CM	
		LH-A201-XL (Extra Large)			A=20CM B=36CM C=34CM	
Foot90°	Ankle joint (90°)	LH-A202-S (Small)			A=16CM B=28CM C=18CM	
		LH-A202-M (Medium)			A=18CM B=33CM C=21CM	
	Mid-to-distal tibia/fibula	LH-A202-L (Large)			A=20CM B=38CM C=24CM	
		LH-A202-XL (Extra Large)			A=21CM B=43CM C=27CM	
Leg	Proximal tibia/fibula	LH-A203-S (Small)			A=39CM B=30CM C=25CM	
		LH-A203-M (Medium)			A=49CM B=34CM C=27CM	
	Distal femur	LH-A203-L (Large)			A=58CM B=42CM C=34CM	
		LH-A203-XL (Extra Large)			A=68CM B=50CM C=41CM	